Instruction for Use Ambu[®] aView[™]









For use by trained clinicians/physicians only.

For in-hospital use.

For use with Ambu[®] visualization devices.

Further instruction details are available online: www.ambu.com

Pat. Pending

Ambu is a registered trademark and aScope and aView are trademarks of Ambu A/S.

Contents	Page
English (Instruction for Use).....	4-17
Български (Указания за ползване).....	18-33
Česky (Návod k použití).....	34-47
Dansk (Brugsanvisning).....	48-61
Deutsch (Bedienungsanleitung).....	62-75
Ελληνικά (δηγίες Χρήσεως).....	76-91
Español (Manual de instrucciones).....	92-105
Eesti (Kasutusjuhised).....	106-118
Suomi (Käyttöohje).....	119-132
Français (Mode d'emploi).....	133-146
Hrvatski (Upute za uporabu).....	147-159
Magyar (Használati útmutató).....	160-173
Italiano (Manuale d'uso).....	174-188
日本語 (使用法).....	189-201
Lietuviškai (Naudojimo instrukcijos).....	202-214
Latviski (Lietošanas instrukcija).....	215-228
Nederlands (Gebruiksaanwijzing).....	229-243
Norsk (Brukerveiledning).....	244-257
Polski (Instrukcja obsługi).....	258-272
Português (Manual de instruções).....	273-287
Português–Brazil (Manual de instruções).....	288-301
Română (Instrucțiuni de utilizare).....	302-316
Русский (Directions for use).....	317-331
Slovenčina (Návod na použitie).....	332-345
Slovenščina (Navodila za uporabo).....	346-360
Svenska (Instruktionshandbok).....	361-374
Türkçe (Kullanım talimatları).....	375-388
中文 (使用说明).....	389-401

1. Important Information – Read Before Use

Read these safety instructions carefully before using the aView monitor. The Instruction for Use may be updated without further notice. Copies of the current version are available upon request. Please be aware that these instructions do not explain or discuss clinical procedures. They describe only the basic operation and precautions related to the operation of the aView monitor.

1.1. Intended Use

The aView™ monitor is a non-sterile, reusable digital monitor, intended to display live imaging data from Ambu visualisation devices.

1.2. Warnings and Cautions

Failure to observe these warnings and cautions may result in patient injury or damage to the equipment. **Ambu is not responsible for any damage to the system or patient injury resulting from incorrect use.**

WARNINGS


1. aView images must not be used as an independent diagnostic of any pathology. Physicians must interpret and substantiate any finding by other means and in the light of the patient's clinical characteristics.
2. Do not use the aView monitor if it is damaged in any way or if any part of the functional check fails.
3. The aView is not to be used when delivering highly flammable anaesthetic gases to the patient. This could potentially cause patient injury.
4. The aView is not to be used in a MRI environment.
5. Do not use aView during defibrillation.
6. When handling the patient do not simultaneously touch the aView power socket or docking connector.
7. Always watch the live endoscopic image on the aView when advancing or withdrawing a visualization device.
8. To avoid risk of electric shock, this equipment must only be connected to a supply mains with protective earth. To disconnect aView from mains remove the mains plug from the wall outlet.
9. Clean and disinfect the aView monitor after each use according to the instructions in section 5. Disconnect aView from any mains power supply, remove any accessories and make sure the aView is completely turned off before cleaning and disinfection.
10. Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.
11. Electronic equipment and the aView may affect the normal function of each other. If the aView is used adjacent to or stacked with other equipment, observe and verify normal operation of both the aView and the other electronic equipment prior to using it. It may be necessary to adopt procedures for mitigation, such as reorientation or relocation of the equipment or shielding of the room in which it is used. Consult the tables in appendix 1 for guidance in placing aView.
12. Be careful to check whether the image on the screen is a live image or a recorded image and verify that the orientation of the image is as expected.
13. Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the system, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

CAUTIONS

1. Have a suitable backup system readily available for immediate use so the procedure can be continued if a malfunction should occur.
2. US federal law restricts this device for sale only by, or on the order of, a physician.
3. Keep aView dry during preparation, use and storage.
4. Pay attention to the battery symbol indicator on the aView monitor. Recharge aView when the battery level is low (see section 4.1). It is recommended that aView is recharged before every procedure and that a charger be readily available during use.
5. Place or hang the aView monitor on a stable support while in use. Dropping aView could damage it.
6. Position the power cord where it is unlikely to be stepped on. Do not place any objects on the power cord.
7. aView is not intended to be repaired. If defect aView shall be discarded.
8. The batteries in aView are not changeable and must only be removed upon disposal.


2. System Description

The aView monitor can be connected with a range of Ambu visualisation devices (cf. section 2.2 for compatible devices) to display the video image from an Ambu visualisation device. The aView monitor is reusable. No modification of this equipment is allowed.

Ambu® aView™	Part numbers:	Bracket (e.g. for attaching the aView to an I.V. pole)	Part numbers:
	405002000 Model no. JANUS2-W08-R10 (SW versions v2.XX)		405017700

For aView model no., please check backside label on aView.

aView is not available in all countries. Please contact your local sales office.

Power supplies		Part numbers:
 x5	aView power supply manufacturer: FSP Group Inc. aView power supply part number: FSP030-REAM	EU/CH 405014700X (not DK & UK) UK 405013700X DK 405012700X AUS & NZ 405015700X US 405016700X

Accessories		Part numbers:
	aView adaptor cable for composite connection	405000712





2.2. Compatible Devices

The compatible Ambu visualisation devices for aView are:

Green connection port (see "5a" on page 2) :	Blue connection port (see "5b" on page 2):
- Ambu® aScope™ 2	- Ambu® aScope™ 3 Family
	- Ambu® aScope™ 4 Broncho Family

Only one visualization device can be connected to aView at a time. aScope 2, aScope 3 Family and aScope 4 Broncho Family are not available in all countries. Please contact your local sales office.


2.3. aView Parts

No.	Part	Function	Material
1	Monitor Casing	-	PC / ABS Rubber
2	Touch Screen	Displays the image from the Ambu visualization device and a touch screen interface.	PET / Glass
3	Stand	To place the display on a solid surface.	ZDC3
4	Connection for Ambu® visualization devices	Power supply and data connection Protected by a rubber cover.	PC / ABS Rubber
5	Power	Power inlet for charging aView. Protected by a rubber cover.	Rubber
6	Input/output connections	JACK connector (video-out) and USB port.	-
7	Power button	Push button for power ON before procedure and power OFF after procedure.	Rubber
8	Docking connector	For future use.	-
9	Power supply	Powers the system Power cord with country-specific plug.	PC PVC
10	Hexagonal key	To tighten the bolt on the bracket.	CrV
11	Bracket	Secures the monitor to e.g. an IV pole.	POM
12	Pouch hook	Slide the hook through the 3 holes on the bracket. It can now be used to hold visualization device pouches with a hole in top corner for easy storage before and during use.	PAA
13	Adapter cable	For the composite video-out interface (see no. 6).	PE/PVC

Abbreviations: PC (Polycarbonate), PE (Polyethylene), ABS (Acrylonitrile Butadiene Styrene), ZDC3 (Zinc alloy), PVC (Polyvinylchloride), POM (Polyoxymethylene), PAA (Polyarylamide), CrV (Chrome Vanadium)

3. Explanation of Symbols Used

Symbols for the aView	Indication	Symbols for the aView	Indication
	Consult Instruction for Use.		Date of Manufacture, followed by YYYY-MM-DD.
	Caution.		Manufacturer.
	Humidity limitation: relative humidity between 30 and 85% in operating environment.		CE mark. The product complies with the EU Council directive concerning Medical Devices 93/42/EEC
	Atmospheric pressure limitation: between 80 and 109 kPa in operating environment		Waste Bin symbol, indicating that waste must be collected according to local regulation and collection schemes for disposal of batteries.
	Connection to external monitor.		Waste Bin symbol, indicating that waste must be collected according to local regulation and collection schemes for disposal of electronic and electrical waste (WEEE).
	Direct current.	Li-ion	Battery type Lithium ion. Only applicable for the battery inside the aView.
	Alternating current.		Re-chargeable battery. Only applicable for the battery inside the aView.
IP30	Protection against solid objects.		Only for indoor use.
	Lot Number, Batch Code.		Tested to comply with FCC Standards - Medical Equipment.
	Reference Number.		

Medical — general medical equipment as to electrical shock, fire and mechanical hazards only in accordance with
 ANSI/AAMI ES60601-1:2005/(R):2012
 CAN/CSA-C22.2 No. 60601-1:08 (R2013)
 IEC 60601-2-18:2009. UL60601-1, 1st ed.: 2006
 CAN/CSA-C22.2 No. 601.1-M90, 2005
 IEC 60601-2-18:1996. Control no 4UD1.

4. Use of aView

Numbers in gray circles below refer to illustrations on page 2.

4.1. Preparation and Inspection of aView

1. Closely examine aView and all parts for any damage (free from wear and tear) **1**.
2. Place aView on a solid flat surface by using the stand on the back of aView **2a**.
If needed, aView can be placed on a pole by using the supplied bracket **2b**.
3. Switch aView on by pressing the power button **3**. A white hourglass indicating that aView is loading the User Interface will appear on the screen.
4. Check the battery indicator on aView. A full charge lasts approximately 3 hours. Charge the aView monitor if necessary by connecting the aView power supply to the wall socket and insert the power plug into the power inlet of aView **4**. Make sure the power supply is present and working at any time. It is recommended to locate the nearest wall socket before start of the procedure.

The battery icon remains white until one block is left, after which it turns red.
 When remaining battery capacity is 10% the red battery icon starts flashing.

Max. battery status of the aView

Min. battery status of the aView

Fully charged battery still connected to charger

Charging is shown with blocks flashing. Current capacity is shown with non-flashing blocks.

Battery is charging

Battery current capacity

Battery damaged

The power button will light up green when aView is turned ON and not connected to mains and orange when aView is connected to mains. When aView is turned off, fully charged and connected to mains the button will turn green.

If the battery is critically low, the recording functions will become unavailable.



If the aView is turned off and the battery is critically low, the power button will flash orange five times every ten seconds to remind that charging is needed.

5. Connect the selected Ambu visualization device to aView by plugging the cable connector into the corresponding female connector on aView (push aside the rubber cover) **5a** or **5b**. Align the arrows on the visualisation device and aView before insertion.
6. Verify that a live video image appears on the screen. Point the distal end of the Ambu visualization device towards an object, e.g. the palm of your hand **6**.
7. Adjust the image preferences on aView if necessary – refer to the online aView User Guide on www.ambu.com.

4.2. Installing, Preparing and Operating the Ambu Visualization Device

Please refer to the Instruction for Use for the specific visualization device.

4.3. Operating the aView Monitor

Startup Image Mode	User Interface
<p>Starts immediately after the power button is pressed and continues until the User Interface is loaded (after approximately 1 min).</p> <p>A live image from the plugged in Ambu visualization device is available while aView is loading. If no device is connected the screen will be blue.</p>  <p>The screen layout may vary depending on the software version.</p>	<p>A live image is available and the User Interface displays advanced user functions (refer to the online aView User Guide on www.ambu.com).</p> <p>aView starts up in the Blue tab for Live Image Viewing and recording live image.</p>  <p>The screen layout may vary depending on the software version.</p>

Interface nomenclature	
 Power button	 Video recording
 Brightness adjustment	 Viewing and recording live image.
 Contrast adjustment	 Managing saved files.
 Snap shot	 System settings and User Accounts.

Live image vs. recorded image

Round BLUE buttons are shown in the blue Live Image tab  and indicates a live image. Square YELLOW or GREEN buttons are shown in the yellow File Management tab  and the green Settings tab  and indicates a recorded image.

Battery status: During start up, aView powers up and configures the visualization device. If the aView battery icon on the screen changes from fully charged to low battery (red battery) within 30 minutes, aView must be replaced.

Connection to external monitor

The image from an Ambu visualization device can be viewed on an external monitor using the video out interface on aView. Connect the external monitor to the composite interface on the right side of aView using the adaptor cable supplied (refer to the online aView User Guide on www.ambu.com). Please consult the manual of the external monitor for further information on how to connect an external video source by composite.

4.4. After Use

Final steps

1. Disconnect the Ambu visualization device from aView. For disposal of the visualisation device, refer to the Instruction for Use of the specific device.
2. Switch off aView by pressing the power button  for at least 2 seconds. A blue hourglass indicating that aView is powering down will appear on the screen and aView will power down.
3. Clean and disinfect aView (cf. section 5).
4. If the battery level of aView is low, charge the aView (cf. section 4.1).

5. Cleaning and Disinfection of aView

The aView should be cleaned and disinfected before and after each use. It is recommended that aView is cleaned and disinfected before and after use following the instructions below. Ambu has validated these instructions in accordance with AAMI TIR 12 & 30. Any deviation from the instructions should be properly evaluated for effectiveness and potential adverse consequences by the reprocessor, to ensure that the device continues to fulfil its intended purpose.

Clean and disinfect the aView according to good medical practice using one of the procedures below:

Procedure 1 - SANI-CLOTH® BLEACH from PDI

Cleaning

Use a wipe to remove heavy soil. All blood and other body fluids must be thoroughly cleaned from surfaces and objects before disinfection by germicidal wipe.

Disinfecting

1. For heavily soiled surfaces, use a wipe to pre-clean aView prior to disinfecting.
2. Unfold a clean wipe and thoroughly wet the surface of aView.
3. The treated surfaces must remain visibly wet for a full four (4) minutes. Use additional wipes if needed to assure continuous 4 minutes wet contact time.
4. Let aView air dry.

Procedure 2 - SUPER SANI-CLOTH® from PDI

Cleaning

Use a wipe to remove heavy soil. All blood and other body fluids must be thoroughly cleaned from surfaces and objects before disinfection by germicidal wipe.

Disinfecting

1. For heavily soiled surfaces, use a wipe to pre-clean aView prior to disinfecting.
2. Unfold a clean wipe and thoroughly wet the surface of aView.
3. The treated surfaces must remain visibly wet for a full two (2) minutes. Use additional wipes if needed to assure continuous 2 minutes wet contact time.
4. Let aView air dry.

Procedure 3

Cleaning

1. Prepare a cleaning solution using a standard enzymatic detergent prepared per manufacturers recommendations. Recommended detergent: enzymatic, mild pH: 7-9, low foaming (Enzol or equivalent).
2. Soak a sterile gauze in the enzymatic solution and make sure that the gauze is moist and not dripping.
3. Thoroughly clean the button, rubber covers, screen, external casing of the monitor and stand with the moist gauze. Avoid getting the device wet to prevent damaging internal electronic components.
4. Using a sterile soft bristled brush that has been dipped in the enzymatic solution, brush the button until all evidence of soil is removed.
5. Wait for 10 minutes (or the time recommended by the manufacturer of the detergent) to allow the enzymes to activate.
6. Wipe the aView clean using sterile gauze that has been moistened with RO/DI water. Ensure all traces of the detergent are removed.
7. Repeat steps 1 to 6.

Disinfecting

Wipe the surfaces of aView for approximately 15 minutes using a piece of sterile gauze moistened with the alcohol mixture indicated below (approximately once every 2 minutes). Follow safety procedures for the handling of isopropyl. The gauze should be moist and not dripping since liquid can affect the electronics inside the aView. Pay close attention to the button, rubber covers, screen, external casing and stand, slots and gaps on the aView. Use a sterile cotton swab for these areas.

Solution: Isopropyl (alcohol) 95%; Concentration: 70-80%; Preparation: 80cc of 95% Isopropyl (alcohol) added to 20cc of purified water (PURW) (Alternatively, use EPA-registered hospital disinfection wipes containing at least 70% isopropyl. Safety precautions and directions of use of the manufacturer must be followed).

After cleaning and disinfection, the aView must be submitted to the pre-check procedure in section 4.1.

Between procedures, aView must be stored in accordance with local guidelines.

6. Maintenance and Disposal

6.1. Maintenance of Battery

To prolong battery life it is recommended to fully charge the monitor at least every third month and store it in a cool place. If the battery is flat the procedure can take up to 5 hours. The battery should be charged at temperatures between 10 - 40°C.

6.2. Disposal

At the end of product life open up the aView and dispose of the batteries and aView separately in accordance with local guidelines.

7. Technical Product Specifications

7.1. Standards Applied

The aView function conforms with:

- Council Directive 93/42/EEC concerning Medical Devices.
- IEC 60601-1 ed 2 Medical electrical equipment – Part 1: General requirements for safety.
- IEC 60601-1 ed 3.1 Medical electrical equipment - Part 1: General requirements for basic safety and essential performance.
- IEC 60601-1-2: Medical electrical equipment – Part 1-2 General requirements for safety – Collateral standard: Electromagnetic compatibility - Requirements for test.

The Ambu aView power supply conforms with:

- Council Directive 93/42/EEC concerning Medical Devices.
- IEC 60601-1 ed 2 Medical electrical equipment – Part 1: General requirements for safety.
- IEC 60601-1 ed 3.1 Medical electrical equipment - Part 1: General requirements for basic safety and essential performance.

7.2. aView Specifications

Display

Max. resolution	800 * 480
Orientation	Landscape
Display type	8.5" colour TFT LCD
Brightness control	Yes, ("+" / "-")
Contrast control	Yes, ("+" / "-")
Start up time	Immediately after pressing the power button.

Connections

USB connection	Type A
Analog composite video out	RCA connection (use adapter cable included)

Storage and transportation

Temperature	10 ~ 40°C (50 ~ 104°F)
Relative humidity	30 ~ 85%
Atmospheric pressure	80-109 kPa

Memory

Storage capacity	8GB
------------------	-----

Mounting interface

Mounting interface standard	VESA MIS-D, 75 C, VESA FDMI compliant display, Part D, with centre located mounting interface
-----------------------------	---

Electrical power

Power requirement	18V 1.67A DC input
Battery type	10.8V 4300mAh
Battery Operation	Typical battery runtime of a new, fully charged battery (aView turned on and scope connected) is min. 3 hours
Protection against electrical shock	Internally Powered

Fixture

Mounting interface	75 mm (2.96")
Fits poles with thicknesses	10 mm ~ 45 mm (0.4 ~ 1.8")

Operating environment

Temperature	10 ~ 40°C (50 ~ 104°F)
Relative humidity	30 ~ 85%
IP Protection Classification System	The aView is classified IP30.
Atmospheric pressure	80-109 kPa
Altitude	≤ 2000 m

Dimensions

Width	241 mm (9.49")
Height	175 mm (6.89")
Thickness	33.5 mm (1.32")
Weight	1500 g (331lbs)

7.3. aView Power Supply

Dimensions

Weight 360g (0.79 lbs)

Electrical power

Power requirement 100 - 240V AC; 50-60Hz; 0.6A

Power out 18V DC; 1.67A

Protection against electrical shock Class I

Operating environment

Temperature 10 ~ 40° C (50 ~ 104° F)

Storage

Temperature 10 ~ 40°C (50 ~ 104°F)

Relative humidity 10 ~ 90%

Plugs

Between the power supply and aView Ø5.5mm DC jack connector

6 interchangeable types

1. Model NEMA 5 AC grounded power plug
2. Australian configuration: AS3112, AC grounded power plug
3. UK configuration: BS1363, AC grounded power plug
4. European configuration: CEE 7, AC grounded power plug
5. Danish configuration: 2-5a, AC grounded power plug
6. Swiss configuration: Type J, AC grounded power plug

Contact Ambu for further information.

8. Trouble Shooting

If problems occur with the aView monitor, please use this trouble-shooting guide to identify the cause and correct the error.

Problem	Possible cause	Recommended action
No live image on the screen but User Interface is present on the display or the image shown is frozen.	An Ambu visualization device is not connected to aView.	Connect an Ambu visualization device to the blue or green port on aView.
	aView and the Ambu visualization device have communication problems.	Restart aView by pressing the power button for at least 2 seconds. When aView is off restart by pressing the power button once more.
	The Ambu visualization device is damaged.	Replace the Ambu visualization device with a new one.
	A recorded image is shown in the yellow file management tab.	Return to live image by pressing the blue Live image tab or restart aView by pressing the power button for at least 2 seconds. When aView is off restart by pressing power button once more.
Low picture quality.	Light reflecting on the aView screen.	Move aView to a position where no direct light influences the screen.
	Dirty / damp screen.	Wipe the screen with a clean cloth.
	Brightness and contrast settings not optimal.	Adjust the contrast and brightness using the designated menu on aView.
Video out is not working.	No live image on the external screen.	Make sure the external monitor is connected to aView using the composite cable and that the external monitor is displaying the correct input. The external monitor shall be capable of receiving NTSC or PAL. Toggle between NTSC and PAL in the video out menu on aView for best result. The external monitor is only displaying a live image when aView is displaying a live image. The image size on the external monitor is determined by the settings in the external monitor.
	The colours, brightness and contrast looks different from the aView screen.	Adjust the colours, brightness and contrast on the external monitor to achieve the desired result.
	Horizontal lines (noise) on external monitor.	If horizontal lines (noise) appear on the external monitor, try to connect the aView power supply.

Appendix 1. Electromagnetic Compatibility

Like other electrical medical equipment the system requires special precautions to ensure electromagnetic compatibility with other electrical medical devices. To ensure electromagnetic compatibility (EMC) the system must be installed and operated according to the EMC information provided in this manual.

The system has been designed and tested to comply with IEC 60601-1-2 requirements for EMC with other devices.

Guidance and manufacturer's declaration – electromagnetic immunity

The system is intended for use in the electromagnetic environment specified below. The customer or the user of the system should assure that it is used in such an environment.

Emissions Test	Compliance	Electromagnetic Environment Guidance
RF emission CISPR 11	Group 1	The system use RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emission CISPR 11	Class B	The system is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emission IEC/EN 61000-3-2	Not applicable	
Voltage fluctuations / flicker emissions IEC/EN 61000-3-3	Complies	

Guidance and manufacturer's declaration – electromagnetic immunity

The system is intended for use in the electromagnetic environment specified below. The customer or the user of the system should assure that it is used in such an environment.

Immunity Test	IEC 60601-1 test level	Compliance Level	Electromagnetic Environment Guidance
Electrostatic discharge (ESD) IEC 61000-4-2	+/- 8 kV contact +/- 2, 4, 8, 15 kV air	+/- 8 kV contact +/- 2, 4, 8, 15 kV air	If floors are covered with synthetic material the relative humidity shall be least 30 %.
Electrical fast transient / burst IEC 61000-4-4	+/- 2 kV for power supply lines +/- 1 kV for input / output lines	+/- 2 kV power supply lines N/A	Mains power quality shall be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	+/- 1 kV line(s) to line(s) +/- 2 kV line(s) to earth		Mains power quality shall be that of a typical commercial or hospital environment.
Voltage Dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5% Ut (95% dip in Ut) for 0.5 cycle 40% Ut (60% dip in Ut) for 5 cycles 70% Ut (30% dip in Ut) for 25 cycles <5% Ut (95% dip in Ut) for 5 sec.	100% reduction 0.5 period 40% reduction for 5 periods 30% reduction for 25 periods 100% reduction for 5 sec.	Mains power quality shall be that of a typical commercial or hospital environment. If the use of the system requires continued operation during power mains interruptions the system can be powered by the built in rechargeable battery.

Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
--	--------	--------	---

Guidance and manufacturer's declaration – electromagnetic immunity

The system is intended for use in the electromagnetic environment specified below.
The customer or the user of the system should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic Environment Guidance
Conducted Radio Frequency IEC 61000-4-6	3 V RMS 0,15 MHz – 80 MHz 6 V RMS in ISM bands 80 % AM at 1 kHz	3 V RMS 0,15 MHz – 80 MHz 6 V RMS in ISM bands 80 % AM at 1 kHz	Portable and mobile RF communications equipment should be used no closer to any part of the system, including its cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.
Radiated Radio Frequency IEC 61000-4-3	3 V/m 80 MHz to 2.7 GHz 80 % AM at 1 kHz	3 V/m 80-2700 MHz 80% AM at 1 kHz	Recommended separation distance $d = 1.17\sqrt{P}$ $d = 1.17\sqrt{P}$ 80 MHz to 800 MHz $d = 2.33\sqrt{P}$ 800 MHz to 2.7 GHz Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, a) should be less than the compliance level in each frequency range b). Interference may occur in the vicinity of equipment marked with the following symbol. 

NOTE 1: At 80 MHz, the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

a) Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast, and TV broadcast, cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey shall be considered. If the measured field strength in the location in which the system is used exceeds the applicable RF compliance level above, the system shall be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the system.

b) Over the frequency range 150kHz to 80MHz, field strengths shall be less than 3 V/m.

Recommended Separation Distances Between Portable and Mobile RF Communication Equipment and system.			
The system is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The user of the system can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters and the system as recommended below, according to the maximum output power of the communication equipment.			
Rated maximum output power (W) of transmitter	Separation distance (m) according to frequency of transmitter		
	150kHz to 80MHz $D = 1.17\sqrt{P}$	80MHz to 800MHz $D = 1.17\sqrt{P}$	800MHz to 2.7GHz $D = 2.33\sqrt{P}$
0.01	0.12 m	0.12 m	0.23 m
0.1	0.37 m	0.37 m	0.74 m
1	1.17 m	1.17 m	2.33 m
10	3.70 m	3.70 m	7.37 m
100	11.7 m	11.7 m	23.3 m

For transmitters rated at a maximum output power not listed above, the recommended separation distance (D) in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Appendix 2. Warranty and Replacement Program

The warranty period for the aView is one year from delivery to the customer. Ambu agrees to replace an aView free of charge if proof can be provided of faulty materials or faulty workmanship. In doing so Ambu cannot accept the cost of transportation or risk of shipment. Repair of an aView due to faulty materials or faulty workmanship may be offered in situations where this option is available. Ambu holds the right to make the choice between repair or replacement in every situation.

Ambu holds the right to ask the technical department at the customer site, or a similar department, to carry out repair of an aView under proper guidance by Ambu. A defective aView must be handled exclusively by persons authorised by Ambu A/S. To prevent infection, it is strictly forbidden to ship contaminated medical devices. The medical device (aView or Ambu visualization device) must be decontaminated on site before shipment to Ambu. The cleaning and disinfection procedures explained in section 5 shall be followed. Ambu reserve the right to return contaminated medical devices to the sender.

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