

Cricothyrotomy Simulator

DATASHEET

The Cricothyrotomy Simulator is the ideal aid for the study and practice of the technique needed to perform surgical or needle cricothyrotomy procedures - perfect for paramedics, EMTs, combat medics, flight nurses, anesthesiologists and other emergency medical personnel.

The simulator has been designed with anatomically accurate landmarks to aid in site training and give confidence in fast action. The neck is hyperextended to enable the user to determine the incision site.

The trachea can be replaced and the airway runs from top to bottom of the simulator, which enables the stylet and obturator placement after incision.

The model has a chin, full-size neck, and ties can be used to hold the obturator securely. A simulated lung can be inflated to ensure correct positioning.



KEY BENEFITS

- Training of surgical or needle cricothyrotomy
- Adult and child tracheas
- Anatomically accurate landmarks
- Hyperextended neck
- Replaceable trachea
- Complete with chin and full-size neck
- Simulated lung verifies correct placement

ORDERING SPECIFICATIONS

Item no.	
LF00994	Cricothyrotomy Simulator incl. simulator on stand, 6 replaceable neck skins, 6 adult trachea inserts (4 rigid, 2 soft), 6 child trachea inserts (4 rigid, 2 soft), simulated lung, lubricant, hard carry case
LF00994 AU	Skins Neck Replacement ST/6
LF00994 BU	Cricoid Replacement Adult ST/6
LF00994 CU	Cricoid Replacement Child ST/6
101-132U	Adult Airway Pkg of 4 for Cricothyrotomy Simulator
101-133U	3 Year Old Airway Pkg of 4 for Cricothyrotomy Simulator
101-136U	Deluxe Overlay Skins Pkg of 4 for Cricothyrotomy Simulator
101-137U	Adult Trachea for Cricothyrotomy Simulator
101-138U	3 Year Old Trachea for Cricothyrotomy Simulator

Ambu

■ Ambu A/S
Baltorpbakken 13
2750 Ballerup
Denmark
T +45 72 25 20 00
F +45 72 25 20 50
ambu.com